

## AGITADOR DE TAMICES RO-TAP MODELO RX-812



Test Sieves sold separately.

### RX-812 Coarse Sieve Shaker Model

Designed to be used in coarse testing applications, the affordable RX-812 Coarse Sieve Shaker has the unique capability of using both 8" diameter and 12" diameter test sieves, simply by using the quick-change adapters provided. The RX-812 has a quiet operation and utilizes an effective action that gives consistent, accurate results.

- Designed for use with 8" and 12" diameter test sieves.
- Not recommended for analysis with particles 150 micron (100 mesh and finer) and smaller.

#### SPECIFICATIONS

- Oscillations per Minute: 280 (60Hz); 231 (50Hz).

Recommendations regarding the maximum number of test sieves that can be used with Ro-Tap Test Sieve Shakers:

MODEL	NUMBER OF SIEVES		PART #	DESCRIPTION	PRICE (\$)
	FULL HEIGHT	HALF HEIGHT			
RX-812 (8" or 200mm Dia.)	6 Sieves w/FH Pan	13 Sieves w/HH Pan	RX-812	115v. 60 Hz Current	
RX-812 (12" or 300mm Dia.)	4 Sieves w/FH Pan	6 Sieves w/HH Pan* or 8 Sieves w/HH Pan	RX-812-1	230v. 60 Hz Current	
			RX-812-2	115v. 50 Hz Current	
			RX-812-3	230v. 50 Hz Current	

A top cover is required for testing.

\* Intermediate height.

Prices do not include test sieves.

All models are wired for single phase current only.

Dimensions — Unit: 23" wide x 15" deep x 13" high (59cm x 38cm x 33cm)  
Shipping: 25" wide x 19" deep x 15" high (64cm x 49cm x 38cm)

Weight — Unit: 55 lbs. (25 kg)  
Shipping: 95 lbs. (28 kg)